

Make it Simple & Esthetic

# T9

Zirconia Sintering Furnace



**BAOT**<sup>®</sup>

**BAOT BIOLOGICAL TECHNOLOGY CO., LTD.**

**Address:** Unit 1, Second Floor, No.12 Building, Yujing Industry Zone, 106 Qihao Road,  
Torch Development District, Zhongshan, Guangdong, P.R.China

**Tel:** +86-760-87893825

**Fax:** +86-760-87893827

**Web:** [www.baotdent.com](http://www.baotdent.com)

1550°C Max / 18 Programs / 5-stage Curves

# BAOT T9 Zirconia Sintering Furnace

The BAOT T9 is a **new silicon carbide rod** zirconia sintering furnace that minimizes zirconia contamination with advanced heating technology.

It features an LCD for real-time monitoring of voltage, current, and humidity, ensuring efficient and stable sintering.

With a maximum temperature of **1550°C** and **18** programmable memory slots, each supporting up to **5-stage** heating curves, it meets the diverse sintering needs of various zirconia brands, making it a reliable solution.



## Features



**Silicon carbide heating:**  
Reduces zirconia contamination



**LCD display:**  
Real-time voltage, current, and humidity



**High temp & multi-program:**  
1550°C max, 18 programs, 5-stage curves



# Technical Specification

Dimensions	42x46.5x80 cm
Packaging Dimensions	59x56x93 cm
Net Weight	44KG
Gross Weight	61KG
Rated Voltage	220V+10%
Rated Frequency	50Hz
Rated Power	2500W
Protection Level	IP21
Fuse 1	2X~380V/16A
Fuse 2	~250V/0.5A
Fuse 3	250V/1.0A
Fast Fuse	~500V/32A

## Spare Part



Firing Table



Sinter Bowl



Ceramic Sintered Beads



Sintering Tray



Tongs



Air Switch